External Deionized Water QC Worksheet

LifeSouth Community Blood Centers

Directions: Collect samples in sterile containers using aseptic technique. Complete the Location and Time Sampled columns below. Samples must be transported at 2 to 8°C and must be received within 24 hours of collection.	Date Sampled:	Sampled by:
	Send results to (email or fax):	
	Email request form to QC@lifesouth.org prior to sending samples.	

Completed by QC Laboratory

Samples received within 24 hours at 2 to 8°C? Yes No Notify requesting facility if samples are unacceptable.		Date/Time Received:	Tech:
Incubator Temperature: °C	In range (22.5 ± 2.5°C)?	Incubation Date/Time:	Tech:

Comple Submittin		Completed by QC Laboratory			
		24-Hour Reading	48-Hour Reading	72-Hour Reading	
	Time	Reading Tech:	Reading Tech:	Reading Tech:	
Location	Sampled	Date:	Date:	Date:	Interpretation
		Time:	Time:	Time:	
		No Growth	No Growth	No Growth	Satisfactory
		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory
		No Growth	No Growth	No Growth	Satisfactory
□ N/A		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory
		No Growth	No Growth	No Growth	Satisfactory
□ N/A		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory
		No Growth	No Growth	No Growth	Satisfactory
□ N/A		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory
		No Growth	No Growth	No Growth	Satisfactory
□ N/A		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory
		No Growth	No Growth	No Growth	Satisfactory
□ N/A		Growth: CFU	Growth: CFU	Growth: CFU	Unsatisfactory

Comments/Corrective Action:

Acceptability Criteria

	Satisfactory	Unsatisfactory
DI Water Cultures	≤ 10 CFU	> 10 CFU

Reviewed and Results sent by: _____ Date: _____